


Transverse divider 1/2 StU 575×280 mm

P/N: 7500218 | QT 278/118/008

HUPFER
we make work flow



Technical data

Weight:	0.2 kg
Width:	278 mm
Depth:	8 mm
Height:	118 mm

Similar to illustration, technical modifications reserved. Without decoration.

The cross divider is used for the optimized separation and organisation of medical instruments in sterilisation baskets.

The cross divider enables a secure and organised arrangement of your medical instruments. Thanks to the effortless hooking and removal, you can easily adjust the cross dividers.

The use of high-quality, recyclable stainless steel ensures long-lasting use and hygienic cleaning.

- Creates clarity
- Use of high-quality, recyclable stainless steel for long-lasting durability and hygienic cleaning

Time and date of the request:
12.02.2026, 13:04:58

All information / dimensions are approximate, technical changes reserved. © Hupfer