

Sterile goods transport cart for Da Vinci 3ERS

P/N: 7506498 | STW/M 1119/769/1062 Da Vinci

HUPFER
we make work flow



Technical data

Insertion type:	Crosswise insertion
TM65 Midlevel Report	Link to the certificate
Payload:	90
Weight:	96 kg
Width:	1119 mm
Depth:	769 mm
Height:	1062 mm

Similar to illustration, technical modifications reserved. Without decoration.

The container transport cart carries and protects sterile goods baskets Da Vinci loaded with Da Vinci endoscope baskets including soft packaging and the prepared instrument.

The efficient Da Vinci stainless steel container transport cart is used for convenient transport and storage of Da Vinci sterile goods baskets with the appropriate Da Vinci endoscope baskets and their soft packaging. The Hupfer sterile goods cart ensures easy handling and efficiency in any medical working environment.

The enclosed design protects the valuable sterile goods from environmental influences and contamination during transport.

The installed drawer frame sets are the perfect fit for Da Vinci sterile goods baskets from Hupfer. Each Da Vinci sterile goods basket has a convenient place for an endoscope basket, its soft packaging and the prepared instrument. The height-adjustable design of the drawer frame sets allows ergonomic and safe removal of the transported instruments.

Thanks to its compact dimensions and low height, the Da Vinci container transport cart provides a good overview and is easy to manoeuvre.

Equipped with two lockable swivel castors and two fixed castors, it moves smoothly

Time and date of the request:
22.02.2026, 09:59:36

All information / dimensions are approximate, technical changes reserved. © Hupfer

Sterile goods transport cart for Da Vinci 3ERS

P/N: 7506498 | STW/M 1119/769/1062 Da Vinci

HUPFER
we make work flow

and stays on track. The antistatic tyres prevent electrostatic charging and effectively protect highly sensitive technical equipment and personnel from potential sparks.

- Da Vinci drawer frame sets allow ergonomic and safe removal of sterile goods
- Enclosed design protects sensitive instruments from environmental influences
- Height-adjustable drawer frame sets allow spacing to be adjusted as necessary.
- Antistatic tyres prevent electrostatic charging and provide effective protection from sparks.
- Swivel castors ensure smooth movement and precise manoeuvring
- Locks guarantee safe movement and immobilisation
- High-quality workmanship in stainless steel enables easy cleaning and perfect hygiene
- Valuable materials ensure sustainability and value retention

Time and date of the request:
22.02.2026, 09:59:36

All information / dimensions are approximate, technical changes reserved. © Hupfer