

## Supply and disposal of sterile goods

**- with the HUPFER® Container und sterile goods transport trolley you close your process chain!**



This transport trolley combines the transport of sterile goods and surgery consumables. Thanks to the re-arranged interior you are even more flexible when transporting instruments or packaging bigger than a StU container or basket.

For even more flexibility our range offers the new CSTWS 9, equipped with two compartments. Solid shelves can be mounted to one of the compartments enabling to store 2 StU side by side or bigger sized sterile goods. The second compartment is equipped with L-shaped insert rails where you can place sterile goods containers or sterile goods baskets. Both the shelves and the rails are adjustable in height inside the trolley so that you can adapt your trolley individually to the goods to be transported.

### Advantages and benefits

- Wide range of equipment and accessories.
- Insert rails and shelves adjustable in height.
- Doors are shockproof in closed and locked open position.
- Vertical push handles for ergonomic handling.



Dimensions W x D x H in mm	Dimensions inside the corpus in mm
1285 x 751 x 1471	1045 x 650 x 1126

**Shelves, L-shaped rails as well as support beams can be removed from the trolley without tools for proper hygienic cleaning.**

All details / dimensions are approximate, we retain the right to make technical changes. © HUPFER® 2012.

HUPFER® Metallwerke GmbH & Co. KG

Dieselstrasse 20 · D-48653 Coesfeld · Postfach 1463 · D-48634 Coesfeld · Tel.: +49 (0) 25 41/805-0 · Fax: +49 (0) 25 41/805-379 · info@hupfer.de · www.hupfer.de