

# One for all: The universal holder for pump bottles

#### Very simple to operate - no hand contact whatsoever

No squeezing half-empty bottles, no dribbling, no hand contact – the new universal holder from Hupfer, for disinfectant and other dispensers, can be conveniently operated with your elbow alone.





### Refilling: clean and fast

Easy and clean refilling is now also guaranteed: flip the holder up, remove the empty bottle, insert a refilled bottle, flip the pump bracket down again, and that's it!



## Flexible height adjustment

The height can be adjusted to the bottle using just two screws. Once you have adjusted the height, it can be fixed and the bottle inserted.

### **Extra flexibility**



Use it for disinfectant, washing lotion or liquid soap: all types of EN bottles, rectangular, round, half-litre and litre bottles, disposable and reusable containers will fit in the holder.

For polymer bottles, e.g.

- rectangular 500 ml (ca. 70 x 55 x 160 mm)
- rectangular 1000 ml (ca. 90 x 75 x 190 mm)
- round 1000 ml
  (ca. Ø 90 x 200 mm)

(measurements always exclude pump height)

## **Technical Data**

Model: Universal hygiene dispenser for bottles with pumping equipment (Art. 8600138)\*

Measurements approx.:120 x 198 x 238-300 mmNet weight:0,88 kgAdjustment range for<br/>dosing lever approx.210 - 270 mm

\*Pump bottle not included!



**Don't give thieves a chance:** As an option, the bottles can be easily secured to the holder using a cable tie.